

Quantitative PCR for detection of *Phytophthora* genus and species-specific: *P. sojae* – *P. sansomeana*

I. Reagents, primers and probes

A. Reagents

- RealMaster Mix Probe without Rox (5 Prime; VWR, Cat No. 10052-276)
 - Includes magnesium acetate 25 mM
- Molecular grade water
- Microcentrifuge amber-colored tubes
- Optical film or optical tubes for qPCR
- Plate or tube strips compatible with the qPCR platform

B. Primers

Primer	Sequence (5' - 3')	Length (bp)	Target
PhyG_ATP9_2FTail	AATAAATCATAACCTTCTTTACAACAAGAATTAATG	36	<i>Phytophthora</i>
PhyG-R6_Tail	AATAAATCATAAATACATAATTCATTTTTATA	32	<i>Phytophthora</i>
FMPI2b	GCGTGGACCTGGAATGACTA	20	Plant IC
FMPI3b	AGGTTGTATTAAAGTTTCGATCG	23	Plant IC

C. Probes

Probes	Sequence (5' - 3')	Length (bp)	Supplier
<i>Phytophthora</i> genus-specific TaqMan probe	[FAM] AAAGCCATC [ZEN] ATTAAACARAATAAAGC [IABkFQ]	26	IDT
<i>P. sojae</i> species-specific TaqMan probe	[HEX] TTGATATAT [ZEN] GAATACAAAGATAGATTTAAGTAAAT [IABkFQ]	35	IDT
<i>P. sansomeana</i> species-specific TaqMan probe	[Quasar670] TATTAGTACTAAYTACTAATATGCATTATTTTTAG [BHQ-2]	35	Biosearch
Plant-IC probe	[CalFluorRed610] CTTTTATTATCACTTCCGGTACTGGCAGG [BHQ-2]	29	Biosearch
Soil-IC (PPF)	[CalFluorRed610] AAAGTAAGCTTATCGATACCGTCGACCT [BHQ-2]	28	Biosearch

* IDT: Integrated DNA Technologies, Inc. (<https://www.idtdna.com/>)

II. Preparation of probes and primers

1. Briefly centrifuge the tubes with primers and probes for 20 secs at 10,000 rpm before opening the tubes.
2. **100 μ M stock solutions:** dissolve primers and probes by adding molecular grade water for a final concentration of 100 μ M. Store probes in amber-colored tubes at -20°C or -80°C if possible.
3. **10 μ M working solutions:** Prepare working solutions by making a ten-fold dilution from the stock solutions using molecular grade water and store 100 μ L aliquots at -20°C or -80°C if possible. Use amber-colored tubes for probes.
 - a. For example, to make 500 μ L of primer working stock, dilute 50 μ L of the primer/probe and dissolve in 450 μ L of molecular grade water.
4. **Plant Internal Control 1 μ M working stocks:** the plant internal control for amplification is required to be at a lower concentration than the rest of the working stocks, hence a 100-fold dilution from the stock solution should be made using molecular grade water. Store 100 μ L aliquots at -20°C or -80°C, using amber-colored tubes for probes.
 - a. For example, to make 500 μ L of plant internal control, dilute 5 μ L of plant internal control primer/probe and dissolve in 495 μ L of molecular grade water.

III. Master mix preparation

1. Thaw primers, probe, RealMaster mix and magnesium acetate. Vortex primers, probe and magnesium acetate (Mg^{++}) briefly for 10 sec at setting 7, then spin for 10 sec in a microspin centrifuge. Place all the reagents on ice, do not vortex the RealMaster mix; mix it by flicking the tube.
2. Prepare the master mix in a clean area; make the mix volume using the table below, adding 5% extra reactions to account for pipetting errors. Keep master mix on ice before use.

A. Plant samples

Reagents	Initial Concentration	Plant samples - Volume per reaction (μL)		Check
		1X	Rxn X	
Primers				
PhyG_ATP9_2FTail	10 μM	1.0		
PhyG-R6_Tail	10 μM	1.0		
Probes				
<i>Phytophthora</i> genus-specific TaqMan probe	10 μM	0.05		
<i>P. sojae</i> species-specific TaqMan probe	10 μM	0.2		
<i>P. sansomeana</i> species-specific TaqMan probe	10 μM	0.1		
Plant Internal Control				
FMPI2b	1 μM	0.4		
FMPI3b	1 μM	0.4		
Plant-IC probe	1 μM	0.4		
Real Master Mix without Rox (5 Prime)	2.5X	8.0		
Mg ⁺⁺	25 mM	2.0		
PCR-grade water		4.45		
Total volume		18 μL		

B. Soil samples

Reagents	Initial Concentration	Soil samples - volume per reaction (μL)		
		1X	Rxn X	Check
Primers				
PhyG_ATP9_2FTail	10 μM	1.0		
PhyG-R6_Tail	10 μM	1.0		
Probes				
<i>Phytophthora</i> genus-specific TaqMan probe	10 μM	0.05		
<i>P. sojae</i> species-specific TaqMan probe	10 μM	0.2		
<i>P. sansomeana</i> species-specific TaqMan probe	10 μM	0.1		
Internal Control (Soil)				
PPF	10 μM	0.2		
Internal Control		1.0		
Real Master Mix without Rox (5 Prime)	2.5X	8.0		
Mg ⁺⁺	25 mM	2.0		
PCR-grade water		4.45		
Total volume		18 μL		

3. Mix the master mix by inverting tube up and down and spinning it down in a microspin centrifuge for 10 sec. Aliquot 18 μ L into each well or tube, keep tubes or plates on ice.
4. Add 2 μ L of sample DNA or 2 μ L of molecular grade water for the NTC control.
 - a. Include 10 pg, 1pg and 100 fg of *Phytophthora sojae* or *Phytophthora sansomeana* as positive controls in each assay.
 - b. The final volume for each reaction after adding DNA should be 20 μ L.
5. Seal the plates or close the tubes and spin down the reaction mix before setting the samples on the platform and keep samples on ice.

IV. qPCR cycling conditions

- Verify the following conditions on the real-time PCR platform:
 1. Stage 1: Denature at 95°C for 2 min
 2. Stage 2: Repeat 45 cycles of:
 - a. 95°C for 15 sec
 - b. 57°C for 90 sec
- Depending on the real-time PCR platform select the appropriate channel configuration:
 - SmartCycler: dye set “FTTC25”
 - StepOne Plus: FAM, VIC, ROX
 - In this case uncheck the ROX as the background fluorescence
 - *P. sansomeana* probe will not work on this platform, so it has to be run in an independent run labeling the probe with HEX.
 - CFX96 or iQ5 systems: FAM, HEX, CalFluor610 and Quasar670

V. Result interpretation

A plant sample should have an internal control with a Ct value in the range of 13 to 30, otherwise samples should be repeated and possibly include a 1:10 diluted sample to improve amplification and help rule out PCR inhibition issues or too high of a concentration of DNA.

Target	Label (Dye)	Expected range
<i>Phytophthora</i> genus	FAM	Positive (Ct < 28)
<i>Phytophthora sojae</i>	HEX	Positive (Ct < 31)
<i>Phytophthora sansomeana</i>	Quasar670	Positive (Ct < 31)
Plant Internal Control	CalFluorRed 610	Ct ranging 13 - 30

VI. References

Bilodeau, G. J., Martin, F. N., Coffey, M. D. & Blomquist, C. L. Development of a multiplex assay for genus- and species-specific detection of *Phytophthora* based on differences in mitochondrial gene order. *Phytopathology* 104, 733–748 (2014).

Bilodeau, G. J., Koike, S. T., Uribe, P. & Martin, F. N. Development of an Assay for Rapid Detection and Quantification of *Verticillium dahliae* in Soil. *Phytopathology* 102, 331–343 (2012).