

Phytophthora qPCR worksheet – Plant Samples

Master mix

Reagents	Initial Concentration	Volume per reaction (μL)		
		1X	18 X	Check
<i>Primers</i>				
PhyG_ATP9_2FTail	10 μM	1.0	18	
PhyG-R6_Tail	10 μM	1.0	18	
<i>Probes</i>				
<i>Phytophthora</i> genus-specific probe	10 μM	0.05	0.9	
<i>P. sojae</i> species-specific probe	10 μM	0.2	3.6	
<i>P. sansomeana</i> species-specific probe	10 μM	0.1	1.8	
<i>Plant Internal Control</i>				
FMPI2b	1 μM	0.4	7.2	
FMPI3b	1 μM	0.4	7.2	
Plant-IC probe	1 μM	0.4	7.2	
Real Master Mix without Rox (5 Prime)	2.5X	8.0	144	
Mg ⁺⁺	25 mM	2.0	36	
PCR-grade water		4.45	80.1	
Sample DNA		2 μL	36 μL	
Total volume		20 μL	360 μL	

Sample set up

16 wells (2 qPCR 8-well strips):

Well	1	Check
A	100pg	
B	100pg	
C	10pg	
D	10pg	
E	1pg	
F	1pg	
G	100fg	
H	100fg	

Well	2	Check
A	Unknown_1	
B	Unknown_1	
C	Unknown_2	
D	Unknown_2	
E	Unknown_3	
F	Unknown_3	
G	NTC	
H	NTC	

SDS – *Fusarium virguliforme* qPCR worksheet

Master mix

Reagent	Working concentration	Final concentration	Volume per reaction (μL)		Check
			1X	18X	
2X Taqman universal master mix	2X	1X	10	180	
F6-3 Primer	20 μM	500nM	0.5	9	
R6 Primer	20 μM	500nM	0.5	9	
Prb3 Probe	10 μM	250nM	0.5	9	
IC F primer(Internal Control)	20 μM	600nM	0.6	10.8	
IC R primer	20 μM	200nM	0.2	3.6	
IC Probe	10 μM	200nM	0.4	7.2	
IC DNA	3k copies/ul	600 copies	0.2	3.6	
BSA	20 mg/ul	200ng/μl	0.2	3.6	
H ₂ O	N/A		4.9	88.2	
Sample DNA	N/A		2	36	
Total Volume			20 μL	360 μL	

Sample set up

16 wells (2 qPCR 8-well strips):

Well	1	Check
A	100pg	
B	100pg	
C	10pg	
D	10pg	
E	1pg	
F	1pg	
G	100fg	
H	100fg	

Well	2	Check
A	Unknown_1	
B	Unknown_1	
C	Unknown_2	
D	Unknown_2	
E	Unknown_3	
F	Unknown_3	
G	NTC	
H	NTC	

Phytophthora genus-specific RPA worksheet

Master mix

Reagents	Initial concentration	RPA <i>Phytophthora</i> genus volume per reaction (μL)		
		1X	4.5X	Check
Genus-specific				
TrnM-F	1 μM	0.50	2.3	
TrnM-R	10 μM	1.45	6.5	
TrnM-Probe	10 μM	0.35	1.6	
<i>Plant Internal Control</i>				
Cox1-IPC-F	10 μM	0.63	2.8	
Cox1-IPC-R	10 μM	0.63	2.8	
Cox1-IPC-Probe	10 μM	0.30	1.4	
RPA Buffer		14.75	66.4	
Water		4.15	18.7	
Crude Plant Extract/DNA		1.00	4.5	
Total Volume*		25 μL	112.5 μL	

* This volume is after adding the magnesium acetate (1.25 μL), but this is added just before the start of the reaction on the cap of each tube.

Sample set up

Add 22.75 μL of master mix on each well:

Well	1	Check
A	<i>Sample 1 – Phyt genus</i>	
B		

Phytophthora sojae species-specific RPA worksheet

Master mix

Reagents	Initial concentration	RPA <i>Phytophthora sojae</i> volume per reaction (μL)		
		1X	4.5X	Check
Species-specific				
Atp9-F	1 μM	0.3	1.1	
Psojae-nad9-R	10 μM	2.1	9.3	
Atp9-Probe	10 μM	0.5	2.0	
RPA Buffer		14.8	66.4	
Water		5.30	23.9	
Crude Plant Extract/DNA		1.0	4.5	
Total Volume*		25 μL	112.5 μL	

* This volume is after adding the magnesium acetate (1.25 μL), but this is added just before the start of the reaction on the cap of each tube.

Sample set up

Add 22.75 μL of master mix on each well:

Well	1	Check
A		
B	<i>Sample 1 – P. sojae</i>	

RPA amplification of *Fusarium virguliforme* with lateral flow detection worksheet

Master mix

Reagents	Initial concentration	1X	8.5 X	___5 X	Check
Rehydration buffer	-	14.75µL	125.375 µL		
Water	-	5.6µL	47.6 µL		
FvF30	10 µM	1.05 µL	8.925 µL		
FvR30	10 µM	1.05 µL	8.925 µL		
FvNfo	10 µM	0.3 µL	2.55 µL		
Crude Plant Extract/DNA		1.00 µL	-		
Total Volume*		25 µL*			

* Total volume is 25 µL after adding the Magnesium acetate, but this is added just before the start of the reaction.

Sample set up

2 wells:

Well	1	Check
A	<i>Sample 1</i>	
B	<i>Sample 2</i>	
C		
D		

Sample set up Bio-Rad (CFX96) platform – *Phytophthora* assay

	Group 4		Group 6				
	1	2	3	4	5	6	...
A	100pg	Unk_1	100pg	Unk_1			
B	100pg	Unk_1	100pg	Unk_1			
C	10pg	Unk_2	10pg	Unk_2			
D	10pg	Unk_2	10pg	Unk_2			
E	1pg	Unk_3	1pg	Unk_3			
F	1pg	Unk_3	1pg	Unk_3			
G	100fg	NTC	100fg	NTC			
H	100fg	NTC	100fg	NTC			

Sample set up ABI (StepOne Plus) platform – SDS assay

	Group 3		Group 5				
	1	2	3	4	5	6	...
A	100pg	Unk_1	100pg	Unk_1			
B	100pg	Unk_1	100pg	Unk_1			
C	10pg	Unk_2	10pg	Unk_2			
D	10pg	Unk_2	10pg	Unk_2			
E	1pg	Unk_3	1pg	Unk_3			
F	1pg	Unk_3	1pg	Unk_3			
G	100fg	NTC	100fg	NTC			
H	100fg	NTC	100fg	NTC			

Sample set up Smart-Cycler (Cepheid) platform – SDS and *Phytophthora* assay

Group 1 - Block A - SDS

	1		2
1	100pg	9	Unk_1
2	100pg	10	Unk_1
3	10pg	11	Unk_2
4	10pg	12	Unk_2
5	1pg	13	Unk_3
6	1pg	14	Unk_3
7	100fg	15	NTC
8	100fg	16	NTC

Group 2 - Block B - Phytophthora

	1		2
1	100pg	9	Unk_1
2	100pg	10	Unk_1
3	10pg	11	Unk_2
4	10pg	12	Unk_2
5	1pg	13	Unk_3
6	1pg	14	Unk_3
7	100fg	15	NTC
8	100fg	16	NTC

Group 7 - Block C - SDS

	1		2
1	100pg	9	Unk_1
2	100pg	10	Unk_1
3	10pg	11	Unk_2
4	10pg	12	Unk_2
5	1pg	13	Unk_3
6	1pg	14	Unk_3
7	100fg	15	NTC
8	100fg	16	NTC