

# Diagnostician training

**Martin Chilvers**  
**Jan Byrne**  
**Janette Jacobs**  
**Alejandro Rojas**

**Jie Wang**  
**Mitch Roth**  
**Tim Miles**

# Workshop outline

## Morning Session - qPCR

8:45 Sample preparation overview

9:15 qPCR overview

9:45 Break

10:00 qPCR setup

11:00 qPCR lectures - SDS and *Phytophthora* specific discussion

12:00 Lunch

## Afternoon Session - Isothermal amplification and RPA

1:00 Review of qPCR results

1:30 RPA overview

2:15 RPA sample prep and reaction setup

3:00 Run RPA reactions

3:30 Break

3:40 Review RPA results


4:00 Q & A session

Workshop wrap up - questions and survey





# Why use molecular techniques?

- Not reliant on morphology
- Very useful if traits are limited or plastic
- Also useful due to convergent evolution
- Useful for resolving relationships
- May support/resolve morphological taxonomy




- *Fusarium* spp.
- Overlapping morphology



# Fundamentals of PCR

<http://www.youtube.com/watch?v=ZmqgRPISg0g>



# PCR






# PCR Thermocycler



# Agarose gel electrophoresis


- DNA is separated by migration through a gel due to voltage gradient, small molecules travel faster than larger ones





# Agarose gel electrophoresis

- Using a DNA ladder of known size we can determine the size of our PCR product



# Real-time PCR instrument (examples)



**Bio-Rad**




**Cepheid**



**Applied Biosystems**


# Real-time qPCR assay

- **Hydrolysis probe chemistry**
  - R: reporter fluorescent molecule
  - Q: fluorescent quencher molecule






# qPCR Serial Dilution Standard Curve



$\Delta R_n$ : Delta normalized reporter

# Unknown Sample Quantification

Ct (cycle threshold)



# Real-time PCR Advantages

- Quantifiable result
- Larger dynamic range  $\sim 7$  orders of magnitude
- Quicker time to result than end-point PCR
- Low-medium sample throughput
  - i.e. single tube reactions, 96/384 well plates, OpenArray - 3072 wells