





Isothermal amplification
P. sojae and *P. sansomeana*

Isothermal Amplification



- Multiple isothermal methods developed
 - ~ 11 different techniques
- Loop-mediated isothermal amplification (LAMP)
- Recombinase Polymerase Amplification
- End-Point and Real-time detection

Isothermal Amplification


	SDA	RCA	LAMP	HDA	RPA
	Strand Displacement	Rolling Circle	Loop Mediated Isothermal	Helicase dependent	Recombinase Polymerase
Template	ssDNA	Circular ssDNA	ssDNA	dsDNA	dsDNA
Performance	~90 min >10 copies	~40 min >10 copies	~60 min 1-10 copies	~100 min 100 copies	~20 min 1-10 copies
Temperature	37°C	30 – 65°C	60 – 65°C	64°C	37 – 42°C
Primers	2-4	1	4 - 6	2	2
Tolerance to biological components	Yes	No	Yes	Yes	Yes
Disadvantage	Less efficient in long targets	Ligation is complicated	Complicated primer design	Complex buffer optimization	Stringent conditions**

Diagnostic qPCR for *P. sojae* and *P. sansomeana*




Assay Background: RPA isothermal

Phytophthora Genus-Specific



Phytophthora Species-Specific



- Optimal primer size ~ 30 bp

Diagnostic RPA for Genus and Species-specific


Phytophthora Genus RPA




+






P. sojae species RPA




P. sansomeana species RPA



RPA process: *P. sojae* and *P. sansomeana*




RPA process: *P. sojae* and *P. sansomeana*




Field Samples: *P. sojae* – *P. sansomeana*

Phytophthora Genus-Specific



Phytophthora Species-Specific



Field Samples: *P. sojae* – *P. sansomeana*

Fields	Isolation	Phytophthora ELISA	RPA <i>Phytophthora</i> Genus	
MIPS2	+	+	+	(2/3)
MIPS3	+	+	+	(3/3)
MIPS4	-	+	+	(3/3)
MIPS5	+	+	+	(3/3)
MIPS6	+	+	+	(3/3)
MIPS7	+	NC	+	(2/3)
MIPS8	+	NC	+	(1/3)
MIPS9	+	+	+	(3/3)
MIPS11	-	+	+	(2/3)
MIPS12	-	+	+	(1/1)
ONPS1	+	+	+	(3/3)
ONPS2	+	+	+	(1/2)
ONPS3	+	+	+	(1/3)
ONPS4	+	+	+	(2/2)
ONPS5	-	+	-	(0/2)
ONPS6	-	+	+	(2/2)