

Isothermal Recombinase Polymerase Amplification


Mitch Roth
Chilvers Lab
Michigan State University

The RPA Cycle


The RPA Cycle

All steps operate at low constant temperature (optimum 37°C)


a. Recombinase / oligonucleotide primer complexes form and target homologous DNA




b. Strand exchange forms a D-loop




c. Polymerase initiates synthesis



d. Parental strands separate & synthesis continues



e. Two duplexes form



Oligonucleotide primers



SSB



Recombinase



Polymerase

Kit components:


- Recombinase
- Polymerase
- Special nuclease
- dNTPs
- Single strand binding proteins

User added components:

- DNA template / target
- Forward primer
- Reverse primer
- Unique probe
- Buffer
- Water


RPA Components

Kit components




RPA Reaction Steps


1. Added primers, probes, and DNA template are allowed to mix with the kit components




RPA Reaction Steps



RPA Reaction Steps




RPA Reaction Steps




RPA Reaction Steps

6. Polymerase uses free dNTPs to extend 3' primer, and make the second complementary strand




RPA Reaction Steps




RPA Reaction Steps


8. Recombinase scans template for primer complement, forms "bubble" or D-loop.




RPA Reaction Steps



RPA Reaction Steps




RPA Reaction Steps




RPA Reaction Steps


This entire process is dynamic, so multiple products form. In this case, the two main products are shown at the bottom.




Interpreting Results – real time



Exposing the probe to light causes it to emit a different wavelength of light. This emitted light is absorbed by an internal quencher. The 3' (right) end of the probe is also blocked so polymerase cannot extend from this position.




The special nuclease present in RPA kits will cleave the probe at a specific site called a THF site. This cleavage separates the fluorophore from the quencher and block.



With the quencher absent, the emitted light from the fluorophore is detectable in a real-time fashion. Also, the new 3' end is unblocked, and polymerase can extend from this free end.

Interpreting Results – real time


Smart-DART
Example output



The y axis is fluorescence values, while the x axis is time.

If fluorescence is increasing over time, this indicates that the probe is finding its compliment sequence and being cut. This is a positive result meaning your target pathogen is present.


Twista
Example output




If fluorescence is not increasing over time, or increasing very minimally, then the probe is not finding its compliment nor being cut. This is a negative result, meaning your target pathogen is absent.

Interpreting Results – end point


In this end-point detection, probe and primer have switched places, and the probe does not have a cleavage site



Interpreting Results – end point




Interpreting Results – end point



Example Results

(+) Control Line →
(+) Test Line B →
(+) Test Line A →



← (+) Control Line
← (+) Test Line B
← (-) Test Line A